

— ABLYBIO, Help Your Research



# RGMB Rabbit pAb

货号: **AYP19566**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	47kDa
实测分子量	47kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Mouse brain,Rat brain
细胞定位	Cell membrane,GPI-anchor,Lipid-anchor,Membrane raft
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 210-453 of human RGM B (NP_001012779.2).
------	--

## 靶点信息

研究背景	
------	--

基因ID	285704
基因名	RGMB
Swiss	Q6NW40 ( <a href="https://www.uniprot.org/uniprotkb/Q6NW40/entry">https://www.uniprot.org/uniprotkb/Q6NW40/entry</a> )
别名	RGMB,DRAGON,RGMB Rabbit pAb,DRG11-responsive axonal guidance and outgrowth of neurite

## 产品验证



Western blot analysis of RGMB expressed in HeLa, Mouse brain, Rat brain using RGMB Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)