

# ERCC3 Rabbit pAb

货号: **AYP19622**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	89kDa
实测分子量	89kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,Mouse kidney,Rat kidney
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 512-782 of human ERCC3 (NP_000113.1).
序列	CPMSPEFYREYVAIKTKKRILLYTMNPNKFRACQFLIKFHERRNDKIIVFADNVFALKEYAIRLNKPIYIGPTSQGERMQILQNF KHNPKINTIFISKVGDTSFDLPEANVLIQISSHGGSSRRQEAQRLGRVLRAKKGMVAEEYN AFFYSLVSQDTQEMAYSTKRQ RFLVDQGYSFKVIKLAGMEEEDLAFSTKEEQQLLQKVLAAATDLDAEEEVVAGEFGSRSSQASRRFGTMSSMSGADDT VYMEYHSSRSKAPSKHVHPLFKRFRK

## 靶点信息

研究背景	This gene encodes an ATP-dependent DNA helicase that functions in nucleotide excision repair. The encoded protein is a subunit of basal transcription factor 2 (TFIIH) and, therefore, also functions in class II transcription. Mutations in this gene are associated with Xeroderma pigmentosum B, Cockayne's syndrome, a nd trichothiodystrophy. Alternative splicing results in multiple transcript variants.
------	---

基因ID	2071
基因名	ERCC3
Swiss	P19447
别名	ERCC3;BTF2;GTF2H;RAD25;TFIIH;TTD2;XPB

## 产品验证



Western blot analysis of ERCC3 expressed in Jurkat, Mouse kidney, Rat kidney using ERCC3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)