

— ABLYBIO, Help Your Research



SLC1A5 Rabbit pAb

货号: **AYP19632**

产品信息

反应	Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	33kDa/36kDa/56kDa
实测分子量	57kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat lung
细胞定位	Cell membrane,Melanosome,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human SLC1A5 (NP_005619.1).
------	--

靶点信息

研究背景	The SLC1A5 gene encodes a sodium-dependent neutral amino acid transporter that can act as a receptor for RD114/type D retrovirus (Larriba et al., 2001 [PubMed 11781704]).[supplied by OMIM, Jan 2011]
------	--

基因ID	6510
基因名	SLC1A5
Swiss	Q15758 (https://www.uniprot.org/uniprotkb/Q15758/entry)
别名	SLC1A5,AAAT,ASCT2,ATBO,M7V1,M7VS1,R16,RDRC,SLC1A5 Rabbit pAb,Baboon M7 virus receptor,RD114/simian type D retrovirus receptor,Sodium-dependent neutral amino acid transporter type 2,Solute carrier family 1 member 5,RDR

产品验证



Western blot analysis of SLC1A5 expressed in Rat lung using SLC1A5 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)