

— ABLYBIO, Help Your Research



# CLIP1 Rabbit pAb

货号: **AYP19686**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:100 - 1:200
理论分子量	156kDa/160kDa/162kDa
实测分子量	182kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,Mouse heart
细胞定位	Cell projection,Cytoplasm,Cytoplasmic side,Cytoplasmic vesicle membrane,Peripheral membrane protein, cytoskeleton,ruffle
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1228-1427 of human CLIP1 (NP_002947.1).
------	---

## 靶点信息

研究背景	The protein encoded by this gene links endocytic vesicles to microtubules. This gene is highly expressed in Reed-Sternberg cells of Hodgkin disease. Several transcript variants encoding different isoforms have been found for this gene.
基因ID	6249
基因名	CLIP1
Swiss	P30622 ( <a href="https://www.uniprot.org/uniprotkb/P30622/entry">https://www.uniprot.org/uniprotkb/P30622/entry</a> )
别名	CLIP1,CLIP,CLIP-170,CLIP170,CYLN1,RSN,CLIP1 Rabbit pAb,Cytoplasmic linker protein 1,Cytoplasmic linker protein 170 alpha-2,Reed-Sternberg intermediate filament-associated protein,Restin

## 产品验证



Western blot analysis of CLIP1 expressed in U-87MG, Mouse heart using CLIP1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))