

— ABLYBIO, Help Your Research



EXOC4 Rabbit pAb

货号: **AYP19737**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	53kDa/110kDa
实测分子量	122kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,Jurkat,Mouse brain,Mouse thymus,Mouse liver,Rat brain
细胞定位	centrosome,cytoplasm,cytosol,exocyst,growth cone membrane,myelin sheath abaxonal region,plasma membrane,synapse
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human EXOC4 (NP_068579.3).
------	---

靶点信息

研究背景	The protein encoded by this gene is a component of the exocyst complex, a multiple protein complex essential for targeting exocytic vesicles to specific docking sites on the plasma membrane. Though best characterized in yeast, the component proteins and functions of exocyst complex have been demonstrated to be highly conserved in higher eukaryotes. At least eight components of the exocyst complex, including this protein, are found to interact with the actin cytoskeletal remodeling and vesicle transport machinery. The complex is also essential for the biogenesis of epithelial cell surface polarity. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.
基因ID	60412
基因名	EXOC4
Swiss	Q96A65 (https://www.uniprot.org/uniprotkb/Q96A65/entry)
别名	EXOC4,SEC8,SEC8L1,Sec8p,EXOC4 Rabbit pAb,Exocyst complex component Sec8,KIAA1699

产品验证



Western blot analysis of EXOC4 expressed in Mouse brain, Rat brain, A-431 using EXOC4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)