

RFPL3 Rabbit pAb

货号: **AYP19742**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	32kDa/35kDa
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,HeLa
细胞定位	cytoplasm,nucleoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-288 of human RFPL3 (NP_006595.1).
序列	MAALFQEASSCPVCSDYLEKPMSECGCTVCLKCINSLQKEPHGEDLLCCCCSMVSQRNKIRPNRQLERLVSHIKELEPKLK KILQMNPRMRKFQVDMTLDADTANNFLISDDLRSVRSGLITQNRQDLAERFDVSVILGSPRFTCGRHYWEVDVGTSTE WDLGVCRESVHCKGKIQLTTELGFWTVSLRDGSRLSASTVPLTFLLVDRKLQRVGIFLDMGMQNVSFDAESGSHVYTF RSVSAEEPLRPFLAPSIPPNGDQGVLSICPLMNSGTTDAPVRPGEAK

靶点信息

研究背景	(Microbial infection Stimulates the activity of Human Immunodeficiency Virus 1/HIV-1 pre-integration complex.
基因ID	10738

基因名	RFPL3
Swiss	O75679
别名	RFPL3

产品验证



Western blot analysis of RFPL3 expressed in U-87MG, HeLa using RFPL3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn