

— ABLYBIO, Help Your Research



# ZFP28 Rabbit pAb

货号: **AYP19782**

## 产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	38kDa/98kDa
实测分子量	99kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-349 of human ZFP28 (NP_001295369.1).
------	--

## 靶点信息

研究背景	May be involved in transcriptional regulation. May have a role in embryonic development.
------	--

基因ID	140612
基因名	ZFP28
Swiss	Q8NHY6 ( <a href="https://www.uniprot.org/uniprotkb/Q8NHY6/entry">https://www.uniprot.org/uniprotkb/Q8NHY6/entry</a> )
别名	ZFP28,mkr5,ZFP28 Rabbit pAb,Krueppel-like zinc finger factor X6,KIAA1431

## 产品验证



Western blot analysis of ZFP28 expressed in Mouse brain,Rat brain using ZFP28 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)