

ZFP28 Rabbit pAb

货号: **AYP19782**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	38kDa/98kDa
实测分子量	99kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-349 of human ZFP28 (NP_001295369.1).
序列	MRGAASASVREPTPLPGRGAPRTKPRAGRGPTVGTATLALPARGRPRSRNGLASKGQRGAAPTGPGRALPSRDTALP QERNKKLEAVGTGIEPKAMSQGLVTFGDVAVDFSQEEWELNPIQRNLYRKVMLENYRNLASLGLCVSKPDVISSLEQ GK EPWTVKRKMTRAWCPDLKAVWKIKELPLKDFCEGKLSQAVITERLTSYNLEYSLLGEHWDYDALFETQPGLVTIKNLAVD FRQQLHPAQKNFCKNGIWENNSDLGSAGHCVAKPDLVSLLEQEKEPMMVKRELTGSLFSGECRRRAWCGKSTARESGAF TKKAASSLIYSRKLYPGSLGLERKLNLSGSSK

靶点信息

研究背景	May be involved in transcriptional regulation. May have a role in embryonic development.
基因ID	140612

基因名	ZFP28
Swiss	Q8NHY6
别名	ZFP28;mkr5

产品验证



Western blot analysis of ZFP28 expressed in Mouse brain,Rat brain using ZFP28 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn