

— ABLYBIO, Help Your Research



FUCA2 Rabbit pAb

货号: **AYP19786**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	19kDa/54kDa
实测分子量	45kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,Mouse kidney,Mouse heart,Mouse stomach,Rat kidney,Rat heart
细胞定位	Secreted
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 118-467 of human FUCA 2 (NP_114409.2).
------	--

靶点信息

研究背景	This gene encodes a plasma alpha-L-fucosidase, which represents 10-20% of the total cellular fucosidase activity. The protein is a member of the glycosyl hydrolase 29 family, and catalyzes the hydrolysis of the alpha-1,6-linked fucose joined to the reducing-end N-acetylglucosamine of the carbohydrate moieties of glycoproteins. This enzyme is essential for Helicobacter pylori adhesion to human gastric cancer cells.
基因ID	2519
基因名	FUCA2
Swiss	Q9BTY2 (https://www.uniprot.org/uniprotkb/Q9BTY2/entry)
别名	FUCA2,dJ20N2.5,FUCA2 Rabbit pAb,Alpha-L-fucoside fucohydrolase 2

产品验证



Western blot analysis of FUCA2 expressed in Mouse kidney,Rat kidney,LO2 using FUCA2 Rabbit pAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blo cking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)