

# SNRPF Rabbit pAb

货号: **AYP19808**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	9kDa
实测分子量	12kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HT-29,293T,LO2,SGC-7901,Mouse brain,Mouse testis
细胞定位	Cytoplasm,Nucleus,cytosol
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-86 of human SNRPF (NP_003086.1).
序列	MSLPLNPKPFLNGLTGKPMVKLKWGMEYKGYLVSDGYMNMQLANTEEYIDGALSGHLGEVLIRCNVLYIRGVVEEEEE DGEMRE

## 靶点信息

研究背景	Plays a role in pre-mRNA splicing as a core component of the spliceosomal U1, U2, U4 and U5 small nuclear ribonucleoproteins (snRNPs, the building blocks of the spliceosome. Component of both the pre-catalytic spliceosome B complex and activated spliceosome C complexes. Is also a component of the minor U12 spliceosome. As part of the U7 snRNP it is involved in histone 3'-end processing.
基因ID	6636

基因名	SNRPF
Swiss	P62306
别名	SNRPF;SMF;Sm-F;snRNP-F

## 产品验证



Western blot analysis of SNRPF expressed in Mouse brain, HeLa, HT-29 using SNRPF Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)