

— ABLYBIO, Help Your Research



SMARCD2 Rabbit pAb

货号: **AYP19830**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200
理论分子量	55kDa/57kDa/58kDa
实测分子量	60KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse thymus,Rat thymus
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 454-531 of human SMA RCD2 (NP_001091896.1).
------	---

靶点信息

研究背景	The protein encoded by this gene is a member of the SWI/SNF family of proteins, whose members display helicase and ATPase activities and which are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The encoded protein is part of the large ATP-dependent chromatin remodeling complex SNF/SWI and has sequence similarity to the yeast Swp73 protein. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	6603
基因名	SMARCD2
Swiss	Q92925 (https://www.uniprot.org/uniprotkb/Q92925/entry)
别名	SMARCD2,BAF60B,CRACD2,PRO2451,Rsc6p,SGD2,SMARCD2 Rabbit pAb,60 kDa BRG-1/Brm-associated factor subunit B,BRG1-associated factor 60B

产品验证



Western blot analysis of SMARCD2 expressed in Mouse thymus,Rat thymus using SMARCD2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)