

DEFB121 Rabbit pAb

货号: **AYP19835**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC IF/ICC
推荐浓度	IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:100
理论分子量	8kDa
实测分子量	8kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-76 of human DEFB121 (NP_001011878.1).
序列	MKLLLLLLVTLLLAQVTPVMKCVGKSGRCRTTCKESEVYILCKTEAKCCVDPKYVPVKPKLTDNTSLESTSAV

靶点信息

研究背景	This gene encodes a member of the beta subfamily of defensins. Beta-defensins are antimicrobial peptides that protect tissues and organs from infection by a variety of microorganisms. This gene is found in a cluster with other beta-defensin genes on the long arm of chromosome 20. Two transcript variants encoding different isoforms have been found for this gene.
基因ID	245934

基因名	DEFB121
Swiss	Q5J5C9
别名	DEFB121;DEFB21;ESC42RELC

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn