

— ABLYBIO, Help Your Research



# POLR1B Rabbit pAb

货号: **AYP19851**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	104kDa/107kDa/121kDa/128kDa/132kDa
实测分子量	128kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HT-29,293T,Mouse spleen
细胞定位	Nucleus,nucleolus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 670-900 of human POLR1B (NP_061887.2).
------	--

## 靶点信息

研究背景	
------	--

基因ID	84172
基因名	POLR1B
Swiss	Q9H9Y6
别名	POLR1B,RPA135,RPA2,Rpo1-2,POLR1B Rabbit pAb,DNA-directed RNA polymerase I 135 kDa polypeptide

## 产品验证



Western blot analysis of POLR1B expressed in Mouse spleen, HeLa, HT-29 using POLR1B Rabbit pAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)