

— ABLYBIO, Help Your Research



HTRA1 Rabbit pAb

货号: **AYP19919**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	51kDa
实测分子量	51kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse liver,Rat brain
细胞定位	Cytoplasm,Secreted,cytosol
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 22-160 of human HTRA1 (NP_002766.1).
------	--

靶点信息

研究背景	This gene encodes a member of the trypsin family of serine proteases. This protein is a secreted enzyme that is proposed to regulate the availability of insulin-like growth factors (IGFs) by cleaving IGF-binding proteins. It has also been suggested to be a regulator of cell growth. Variations in the promoter region of this gene are the cause of susceptibility to age-related macular degeneration type 7.
基因ID	5654
基因名	HTRA1
Swiss	Q92743 (https://www.uniprot.org/uniprotkb/Q92743/entry)
别名	HTRA1,ARMD7,CADASIL2,CARASIL,HtrA,L56,ORF480,PRSS11,HTRA1 Rabbit pAb,High-temperature requirement A serine peptidase 1,Serine protease 11

产品验证



Western blot analysis of HTRA1 expressed in Mouse brain, Mouse liver, Rat brain using HTRA1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)