

DDX3Y Rabbit pAb

货号: **AYP19923**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:100
理论分子量	32kDa/72kDa/73kDa
实测分子量	73kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,B cells,293T,A-431,HeLa,Mouse testis,Rat testis
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human DDX3Y (NP_001116137.1).
序列	MSHVVVKNDELQQLANLDLNSEKQSGGASTASKGRYIPPLRNREASKGFHDKDSSGWSCSKDKDAYSSFGSRDSR GKPGYFSERGSRSRGRFDDRGRSDYDGIGNRERPGFRFRSGHSRWCDKSVEDDWSKPLPPSERLEQELFSGGNTGI NFEKYDDIPVEATGSNCPPIENFSDIDMGEIIMGNIELTRYTRPTPVQKHAIP

靶点信息

研究背景	The protein encoded by this gene is a member of the DEAD-box RNA helicase family, characterized by nine conserved motifs, included the conserved Asp-Glu-Ala-Asp (DEAD) motif. These motifs are thought to be involved in ATP binding, hydrolysis, RNA binding, and in the formation of intramolecular interactions. This protein shares high similarity to DDX3X, on the X chromosome, but a deletion of this gene is not complemented by DDX3X. Mutations in this gene result in male infertility, a reduction in germ cell numbers, and can result in Sertoli-cell only syndrome. Pseudogenes sharing similarity to both this gene and the DDX3X paralog are found on chromosome 4 and the X chromosome. Alternative splicing results in multiple transcript variants encoding different isoforms.
基因ID	8653
基因名	DDX3Y
Swiss	O15523
别名	DDX3Y;DBY

产品验证



Western blot analysis of DDX3Y expressed in Mouse testis, Rat testis, A-549 using DDX3Y Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn