

HNF4a Rabbit pAb

货号: **AYP19967**

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	51kDa/52kDa
实测分子量	53kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human HNF4a (NP_000448.3).
序列	RQCVVDKDKRNQCRYCRLKKCFRAGMKKEAVQNERDRISTRSSYEDSSLPSINALLQAEVLSRQITSPVSGINGDIRAKKI ASIADVCEMKEQLLVLE

靶点信息

研究背景	The protein encoded by this gene is a nuclear transcription factor which binds DNA as a homodimer. The encoded protein controls the expression of several genes, including hepatocyte nuclear factor 1 alpha, a transcription factor which regulates the expression of several hepatic genes. This gene may play a role in development of the liver, kidney, and intestines. Mutations in this gene have been associated with monogenic autosomal dominant non-insulin-dependent diabetes mellitus type I. Alternative splicing of this gene results in multiple transcript variants encoding several different isoforms.
------	---

基因ID	15378
基因名	
Swiss	
别名	HNF4A;FRTS4;HNF4;HNF4a7;HNF4a8;HNF4a9;HNF4alpha;MODY;MODY1;NR2A1;NR2A21;TCF;TCF14;HNF4a

产品验证



Western blot analysis of HNF4a expressed in Mouse liver using HNF4a Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn