

POLR3A Rabbit pAb

货号: AYP20032

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200
理论分子量	155kDa
实测分子量	156kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human POLR3A (NP_008986.2).
序列	MVKEQFRETDVAKKISHICFGMKSPPEMRQQAHIQVVSKNLYSQDNQHAPLLYGVLDRMGTSKDRPCETCGKNLADCLGHYGYIDLELPCFHVGYFRAVIGILQMICKTCCHIMLSQEKKQFLDYLRPGLTYLQKRGLKISDKCRKKNICHHCAGFNGTVKKCGLLKIIHEKYKTNKKVVDPIVSNFLQSFETA

靶点信息

研究背景	The protein encoded by this gene is the catalytic component of RNA polymerase III, which synthesizes small RNAs. The encoded protein also acts as a sensor to detect foreign DNA and trigger an innate immune response.
基因ID	11128

基因名	POLR3A
Swiss	O14802
别名	POLR3A;ADDH;HLD7;RPC1;RPC155;hRPC155

产品验证



Western blot analysis of POLR3A expressed in 293T using POLR3A Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn