

# NFKBIZ Rabbit pAb

货号: **AYP20079**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:1000 - 1:5000
理论分子量	64kDa/68kDa/78kDa
实测分子量	64KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat lung
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human NFKBIZ (NP_113607.1).
序列	MIVDKLLDDSRGGEGLRDAAGGCGLMTSPLNLSYFYGASPPAAAPGACDASCSVLGSPAGSPGSDSSDFSSASSVSSC GAVESRSRGGARAERQPVEPHMGVGRQQRGPFQGVVRVKNVSKELLLHIRSHKQKASGQAVDDFKTQGVNIEQFRELKN TVSYSGKRKGPDSLSDGPACKRPALLHSQFLTPPTPTPGESMEDVHLNPKQESSADLLQNI

## 靶点信息

研究背景	This gene is a member of the ankyrin-repeat family and is induced by lipopolysaccharide (LPS). The C-terminal portion of the encoded product which contains the ankyrin repeats, shares high sequence similarity with the I kappa B family of proteins. The latter are known to play a role in inflammatory responses to LPS by their interaction with NF-B proteins through ankyrin-repeat domains. Studies in mouse indicate that this gene product is one of the nuclear I kappa B proteins and an activator of IL-6 production. Two transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	64332
基因名	NFKBIZ
Swiss	Q9BYH8
别名	NFKBIZ;IKBZ;INAP;MAIL

## 产品验证



Western blot analysis of NFKBIZ expressed in Rat lung using NFKBIZ Rabbit pAb at 1:1000. Secondary anti body: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% n-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)