

BANK1 Rabbit pAb

货号: AYP20129

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:1000 - 1:2000
理论分子量	73kDa/86kDa/88kDa/89kDa
实测分子量	89kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	B cells,Raji,U-87MG,A-549,22Rv1,HepG2,Mouse liver,Rat liver
细胞定位	perinuclear region of cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 516-785 of human BANK1 (NP_060405.4).
序列	PVANAFQLERPHFTLPGTMVEGQMERSQNWGHPGVRQETGDEPKGEKEKKEEKEQEEEEEDPYTFAEIDDSEYDMILANLSIKKKTGSRFSIINRPPAPTTPRPTSIPPKEETTPYIAQVFQKQTARRQSDDDKFRGLPKKQDRARIESPAFSTLRGCLTDGQEELILLQEKVKNGKMSMDEALEKFKHWQMGKSGLEMIQKEKLRQLRDCIIGKRPEEENVYNKLTIVHHPGGKETAHNENKFINVHFSNKLPARPQVEKEFGCCKKDH

靶点信息

研究背景	The protein encoded by this gene is a B-cell-specific scaffold protein that functions in B-cell receptor-induced calcium mobilization from intracellular stores. This protein can also promote Lyn-mediated tyrosine phosphorylation of inositol 1,4,5-trisphosphate receptors. Polymorphisms in this gene are associated with susceptibility to systemic lupus erythematosus. Alternative splicing results in multiple transcript variants.
------	--

基因ID	55024
基因名	BANK1
Swiss	Q8NDB2
别名	BANK1;BANK

产品验证



Western blot analysis of BANK1 expressed in Mouse liver,Rat liver,B cells using BANK1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn