

— ABLYBIO, Help Your Research



FBXL2 Rabbit pAb

货号: **AYP20134**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	39kDa/47kDa
实测分子量	38-47kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,A-549,HT-1080,SKOV3,Mouse brain,Mouse kidney,Rat brain,Rat kidney
细胞定位	Lipid-anchor,Membrane
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 120-240 of human FBXL 2 (NP_036289.3).
------	--

靶点信息

研究背景	This gene encodes a member of the F-box protein family which is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. The protein encoded by this gene belongs to the Fbls class and, in addition to an F-box, contains 12 tandem leucine-rich repeats. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	25827
基因名	FBXL2
Swiss	Q9UKC9 (https://www.uniprot.org/uniprotkb/Q9UKC9/entry)
别名	FBXL2,FBL2,FBL3,FBXL2 Rabbit pAb,F-box and leucine-rich repeat protein 2,F-box protein FBL2/FBL3

产品验证



Western blot analysis of FBXL2 expressed in Mouse brain,Rat brain,U-87MG using FBXL2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)