

Aurora B Rabbit pAb

货号: **AYP20169**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200
理论分子量	16kDa/34kDa/35kDa/39kDa
实测分子量	40KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,NIH/3T3,C6
细胞定位	Chromosome,Cytoplasm,Midbody,Nucleus,centromere,cytoskeleton,spindle
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Aurora B (NP_004208.2).
序列	MAQKENSYPWPYGRQTAPSGLSTLPQRVLRKEPVTPSALVLMRSNVQPTAAPGQKVMENSSGTPDILTRHFTIDDFEIGRPLGKGFKNVYLAREKKSH

靶点信息

研究背景	This gene encodes a member of the aurora kinase subfamily of serine/threonine kinases. The genes encoding the other two members of this subfamily are located on chromosomes 19 and 20. These kinases participate in the regulation of alignment and segregation of chromosomes during mitosis and meiosis through association with microtubules. A pseudogene of this gene is located on chromosome 8. Alternatively spliced transcript variants have been found for this gene.
------	--

基因ID	9212
基因名	AURKB
Swiss	Q96GD4
别名	AIK2;AIM-1;AIM1;ARK2;AurB;IPL1;PPP1R48;STK12;STK5;aurkb-sv1;aurkb-sv2;AURKB;Aurora B

产品验证



Western blot analysis of Aurora B expressed in HeLa, Jurkat, NIH/3T3 using Aurora B Rabbit pAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn