

— ABLYBIO, Help Your Research



# ARMS2 Rabbit pAb

货号: **AYP20179**

## 产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:200 - 1:1000
理论分子量	11kDa
实测分子量	11kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse eye
细胞定位	Cytoplasm
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-107 of human ARMS2 (NP_001093137.1).
------	--

## 靶点信息

研究背景	This gene encodes a protein that is thought to play a role in diseases in the elderly. Mutations in this gene have been associated with age-related macular degeneration.
------	---

基因ID	387715
基因名	ARMS2
Swiss	P0C7Q2 ( <a href="https://www.uniprot.org/uniprotkb/P0C7Q2/entry">https://www.uniprot.org/uniprotkb/P0C7Q2/entry</a> )
别名	ARMS2,ARMD8,ARMS2 Rabbit pAb

## 产品验证



Western blot analysis of ARMS2 expressed in Mouse eye using ARMS2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)