

— ABLYBIO, Help Your Research



# CEACAM4 Rabbit pAb

货号: **AYP20234**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:1000
理论分子量	25kDa
实测分子量	24kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,A-549
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 36-155 of human CEACAM4 (NP_001808.2).
------	----------------------------------------------------------------------------------------------------------------------

## 靶点信息

研究背景	Granulocyte orphan receptor that acts as an trigger efficient phagocytosis of attached particles.
------	---------------------------------------------------------------------------------------------------

基因ID	1089
基因名	CEACAM4
Swiss	O75871 ( <a href="https://www.uniprot.org/uniprotkb/O75871/entry">https://www.uniprot.org/uniprotkb/O75871/entry</a> )
别名	CEACAM4,CGM7,CGM7_HUMAN,NCA,CEACAM4 Rabbit pAb,Carcinoembryonic antigen CGM7,Carcinoembryonic antigen-related cell adhesion molecule 4,Non-specific cross-reacting antigen W236

## 产品验证



Western blot analysis of CEACAM4 expressed in U-87MG,A-549 using CEACAM4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)