

— ABLYBIO, Help Your Research



# PPIG Rabbit pAb

货号: **AYP20306**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	40kDa/88kDa
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat, HepG2, HeLa
细胞定位	Nucleus matrix, Nucleus speckle
纯化	Affinity purification

## 抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PPIG (NP_004783.2).
------	---

## 靶点信息

研究背景	
------	--

基因ID	9360
基因名	PPIG
Swiss	Q13427 ( <a href="https://www.uniprot.org/uniprotkb/Q13427/entry">https://www.uniprot.org/uniprotkb/Q13427/entry</a> )
别名	PPIG,CARS-Cyp,CYP,SCAF10,SRCP,PPIG Rabbit pAb,CASP10,Clk-associating RS-cyclophilin,Cyclophilin G, Rotamase G

## 产品验证



Western blot analysis of PPIG expressed in Jurkat, HepG2, HeLa using PPIG Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)