

POLR1C Rabbit pAb

货号: **AYP20316**

产品信息

| | |
|-------|--|
| 反应 | Human |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB IHC |
| 推荐浓度 | WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 |
| 理论分子量 | 38kDa/39kDa |
| 实测分子量 | 39kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | K-562,293T,MCF7 |
| 细胞定位 | Nucleus |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|---|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human POLR1C (NP_976035.1). |
| 序列 | MAASQAVEEMRSRVVLGEFGVRNVHTTDFPGNYSYGYDDAWDQDRFEKNFRVDVVHMDENSLEFDMVGIDAAIANAFR RILLAEVPTMAVEKVLVYNNTSIVQDEILAHRLGLIPIHADPRLFEYRNQGDDEEGTEIDTLQFRLQVRCTRNPAAKDSSDPN ELYVNHKVYTRHMTWIPLGNQADLFPEGTIRPVHDDILIA |

靶点信息

| | |
|------|---|
| 研究背景 | The protein encoded by this gene is a subunit of both RNA polymerase I and RNA polymerase III complexes. The encoded protein is part of the Pol core element. Mutations in this gene have been associated with Treacher Collins syndrome (TCS) and hypomyelinating leukodystrophy 11. Alternative splicing results in multiple transcript variants. |
|------|---|

| | |
|-------|---|
| 基因ID | 9533 |
| 基因名 | POLR1C |
| Swiss | O15160 |
| 别名 | POLR1C;AC40;HLD11;RPA39;RPA40;RPA5;RPAC1;RPC40;TCS3 |

产品验证



Western blot analysis of POLR1C expressed in K-562,293T,MCF7 using POLR1C Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn