

— ABLYBIO, Help Your Research



PDLIM7 Rabbit pAb

货号: **AYP20409**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	50kDa
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,Mouse lung
细胞定位	actin cytoskeleton,cytosol,nucleoplasm,Z disc
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 90-190 of human PDLIM 7 (NP_998801.1).
------	----------------------------------------------------------------------------------------------------------------------

靶点信息

研究背景	
------	--

基因ID	9260
基因名	PDLIM7
Swiss	Q9NR12 (https://www.uniprot.org/uniprotkb/Q9NR12/entry)
别名	LMP1,LMP3,PDLIM7,PDLIM7 Rabbit pAb,LIM mineralization protein,Protein enigma,ENIGMA

产品验证



Western blot analysis of PDLIM7 expressed in Jurkat, Mouse lung using PDLIM7 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)