




基因名	LAMP1
Swiss	P11279
别名	LAMP1;CD107a;LAMPA;LGP120

## 产品验证



Western blot analysis of LAMP1 expressed in Mouse lung, Mouse kidney, Rat liver using LAMP1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)