

— ABLYBIO, Help Your Research



CD70 Rabbit pAb

货号: **AYP20527**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IF/ICC
推荐浓度	IF/ICC: 1:50 - 1:200
理论分子量	21kDa/23kDa
实测分子量	21kDa/23kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Membrane,Single-pass type II membrane protein
纯化	Affinity purification

抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 100-193 of human CD70 (NP_001243.1).
------	---

靶点信息

研究背景	The protein encoded by this gene is a cytokine that belongs to the tumor necrosis factor (TNF) ligand family. This cytokine is a ligand for TNFRSF27/CD27. It is a surface antigen on activated, but not on resting, T and B lymphocytes. It induces proliferation of costimulated T cells, enhances the generation of cytolytic T cells, and contributes to T cell activation. This cytokine is also reported to play a role in regulating B-cell activation, cytotoxic function of natural killer cells, and immunoglobulin synthesis.
基因ID	970
基因名	CD70
Swiss	P32970 (https://www.uniprot.org/uniprotkb/P32970/entry)
别名	CD70,CD27-L,CD27L,CD27LG,TNFSF7,TNLG8A,CD70 Rabbit pAb,CD27 ligand,Tumor necrosis factor ligand superfamily member 7

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn//www.ablybio.cn>)