

ETS2 Rabbit pAb

货号: **AYP20537**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	53kDa
实测分子量	53kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-469 of human ETS2 (P15036).
序列	MNDFGIKNDQVAPVANSYRGLTKRQPAFDTFDGLFAVFPNLNEEQTLQEVPTGLDSISHDSANCELPLLTSCAVMS QALKATFSGFKKEQRRLGIPKNPWLWSEQQVCQWLLWATNEFSLVNVNLQRFGMNGQMLCNLGKERFLELAPDFVGD LWEHLEQMIKENQEKTEDQYEENSHLTSVPHWINSNTLGFGEQAPYGMQTNYPKGGLLDSMCPASTPSVLSSEQEFQ MFPKSRLSSVSVTYCSVSQDFPGSNLNLNNSGTPKDHDSPENGADSFESSDLLQSWNSQSSLLDVQRVPSFESFED DCSQLCLNKPTMSFKDYIERSDPVEQKPVIPAAVLGFTGSGPIQLWQFLLELLSDKSCQSFISWTGDGWEFKLADP DEVARRWGKRKNKPKMNYEKLSRGLRYYYDKNIIHKTS GKRYVYRFVCDLQNLGFTPEELHAILGVQPDTE

靶点信息

研究背景	This gene encodes a transcription factor which regulates genes involved in development and apoptosis. The encoded protein is also a protooncogene and shown to be involved in regulation of telomerase. A pseudogene of this gene is located on the X chromosome. Alternative splicing results in multiple transcript variants.
------	---

基因ID	2114
基因名	ETS2
Swiss	P15036
别名	ETS2;ETS2IT1

产品验证



Western blot analysis of ETS2 expressed in human liver cancer, mouse lung, rat testis using ETS2 Rabbit p Ab at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn