

— ABLYBIO, Help Your Research



# NELL2 Rabbit pAb

货号: **AYP20554**

## 产品信息

|       |  |
|-------|--|
| 反应    | Mouse,Rat  |
| 宿主    | Rabbit   |
| 克隆性   | Polyclonal   |
| 预测反应  |  |
| 应用    | WB   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000  |
| 理论分子量 |  |
| 实测分子量 | 91kDa  |
| 形式    | Liquid   |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated   |
| 阳性对照  | Mouse liver,Mouse brain,Mouse kidney,Rat testis,Rat brain  |
| 细胞定位  | cytoplasm,extracellular region   |
| 纯化    | Affinity purification  |

## 抗原信息


|      |  |
|------|--|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 22-300 of human NELL2 (NP_006150.1). |
|------|--|

## 靶点信息

|      |  |
|------|--|
| 研究背景 |  |
|------|--|

|       |  |
|-------|--|
| 基因ID  | 4753   |
| 基因名   | NELL2  |
| Swiss | Q99435 ( <a href="https://www.uniprot.org/uniprotkb/Q99435/entry">https://www.uniprot.org/uniprotkb/Q99435/entry</a> ) |
| 别名    | NRP2,NELL2,NELL2 Rabbit pAb,NEL-like protein 2,Nel-related protein 2   |

## 产品验证



Western blot analysis of NELL2 expressed in Mouse liver,Rat testis,Mouse brain using NELL2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)