

— ABLYBIO, Help Your Research



# CD11a/LFA-1A/ITGAL Rabbit pAb

货号: **AYP20561**

## 产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	119kDa/128kDa/134kDa
实测分子量	180kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jukat,HepG2,HeLa
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-1170 of human CD11a/LFA-1A/CD11a/LFA-1A/ITGAL (P20701).
------	---

## 靶点信息

研究背景	ITGAL encodes the integrin alpha L chain. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. This I-domain containing alpha integrin combines with the beta 2 chain (ITGB2) to form the integrin lymphocyte function-associated antigen-1 (LFA-1), which is expressed on all leukocytes. LFA-1 plays a central role in leukocyte intercellular adhesion through interactions with its ligands, ICAMs 1-3 (intercellular adhesion molecules 1 through 3), and also functions in lymphocyte costimulatory signaling. Two transcript variants encoding different isoforms have been found for this gene.
基因ID	3683
基因名	ITGAL
Swiss	P20701 ( <a href="https://www.uniprot.org/uniprotkb/P20701/entry">https://www.uniprot.org/uniprotkb/P20701/entry</a> )
别名	CD11A,LFA-1,LFA1A,ITGAL,CD11a/LFA-1A/ITGAL Rabbit pAb,CD11 antigen-like family member A,Leukocyte adhesion glycoprotein LFA-1 alpha chain,Leukocyte function-associated molecule 1 alpha chain

## 产品验证



Western blot analysis of CD11a/LFA-1A/ITGAL expressed in Jukat, HepG2, HeLa using CD11a/LFA-1A/ITGAL Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)