

— ABLYBIO, Help Your Research



# IFNA10 Rabbit pAb

货号: **AYP20566**

## 产品信息

反应	Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	
实测分子量	25kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat plasma
细胞定位	extracellular region,extracellular space
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 24-100 of human IFNA10 (NP_002162.1).
------	---

## 靶点信息

研究背景	This gene encodes a protein that belongs to the type I interferon family of proteins, and is located in a cluster of alpha interferon genes on chromosome 9. Interferons are small regulatory molecules that function in cell signaling in response to viruses and other pathogens or tumor cells. This gene is intronless and the encoded protein is secreted.
基因ID	3446
基因名	IFNA10
Swiss	P01566 ( <a href="https://www.uniprot.org/uniprotkb/P01566/entry">https://www.uniprot.org/uniprotkb/P01566/entry</a> )
别名	IFNA10,IFN-alphaC,IFNA10 Rabbit pAb,Interferon alpha-6L,Interferon alpha-C

## 产品验证



Western blot analysis of IFNA10 expressed in Rat plasma using IFNA10 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)