

— ABLYBIO, Help Your Research



# SCGB1A1 Rabbit pAb

货号: **AYP20612**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	9kDa
实测分子量	11kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung
细胞定位	Secreted
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 22-91 of human SCGB1 A1 (NP_003348.1).
------	--

## 靶点信息

研究背景	This gene encodes a member of the secretoglobin family of small secreted proteins. The encoded protein has been implicated in numerous functions including anti-inflammation, inhibition of phospholipase A2 and the sequestering of hydrophobic ligands. Defects in this gene are associated with a susceptibility to asthma.
基因ID	7356
基因名	SCGB1A1
Swiss	P11684 ( <a href="https://www.uniprot.org/uniprotkb/P11684/entry">https://www.uniprot.org/uniprotkb/P11684/entry</a> )
别名	SCGB1A1,CC10,CC16,CCPBP,CCSP,UGB,UP-1,UP1,SCGB1A1 Rabbit pAb,Club cell phospholipid-binding protein,Club cells 10 kDa secretory protein,Secretoglobin family 1A member 1,Urinary protein 1

## 产品验证



Western blot analysis of SCGB1A1 expressed in Mouse lung using SCGB1A1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)