

— ABLYBIO, Help Your Research



# RPP40 Rabbit pAb

货号: **AYP20666**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	
实测分子量	42kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,U-87MG,Mouse testis,Mouse kidney,Mouse liver,Rat lung,Rat thymus
细胞定位	nucleolar ribonuclease P complex,nucleoplasm,nucleus
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 130-210 of human RPP40 (NP_001273061).
------	--

## 靶点信息

研究背景	Component of ribonuclease P, a ribonucleoprotein complex that generates mature tRNA molecules by cleaving their 5'-ends. Also a component of the MRP ribonuclease complex, which cleaves pre-rRNA sequences.
基因ID	10799
基因名	RPP40
Swiss	O75818 ( <a href="https://www.uniprot.org/uniprotkb/O75818/entry">https://www.uniprot.org/uniprotkb/O75818/entry</a> )
别名	RPP40,RNASEP1,bA428J1.3,RPP40 Rabbit pAb,RNase P subunit 1

## 产品验证



Western blot analysis of RPP40 expressed in Mouse testis,Rat lung,U-251MG using RPP40 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))