

# ASB4 Rabbit pAb

货号: **AYP20692**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC
推荐浓度	<b>IHC:</b> 1:50 - 1:200
理论分子量	
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	cytosol
纯化	Affinity purification

## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human ASB4 ( NP_665879.1).
序列	MDGTTAPVTKSGAAKLVKRNFLKSNDFGKLIKALIQRIDVDTVFEVEDENMVLASYKQGYWLPYKLSKSSWATGLHLSVLFVGHVECLLVLLDHNATINCRPNKGPLHVACEMANVDCVKILCDRGAKLNCYLSGHTALHFCTTPSSILCAKQLVWRGANVNMKTNNQDEETPLHTAAHFGLSELVAFYVEHGAIVDSVNAHMETPLAIAAYWALRFKEQEYSTEHHLVCRMLLDYKAEVNARDDDFKSPLHKAAWNCNHVLMHMMLEAGAEANLMDINGCAAQYVLKVTSVRP

## 靶点信息

研究背景	The protein encoded by this gene is a member of the ankyrin repeat and SOCS box-containing (ASB) family of proteins. They contain ankyrin repeat sequence and SOCS box domain. The SOCS box serves to couple suppressor of cytokine signalling (SOCS) proteins and their binding partners with the elongin B and C complex, possibly targeting them for degradation. Multiple alternatively spliced transcript variants have been described for this gene but some of the full length sequences are not known.
------	--

基因ID	51666
基因名	ASB4
Swiss	Q9Y574
别名	ASB4;ASB-4

产品验证

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)