

— ABLYBIO, Help Your Research



# ZSCAN2 Rabbit pAb

货号: **AYP20700**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	
实测分子量	75kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	nucleus
纯化	Affinity purification


## 抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human ZSCAN2 (NP_870992.2).
------	--

## 靶点信息

研究背景	The protein encoded by this gene contains several copies of zinc finger motif, which is commonly found in transcriptional regulatory proteins. Studies in mice show that this gene is expressed during embryonic development, and specifically in the testis in adult mice, suggesting that it may play a role in regulating genes in germ cells. Alternative splicing of this gene results in several transcript variants encoding different isoforms.
基因ID	54993
基因名	ZSCAN2
Swiss	Q7Z7L9 ( <a href="https://www.uniprot.org/uniprotkb/Q7Z7L9/entry">https://www.uniprot.org/uniprotkb/Q7Z7L9/entry</a> )
别名	ZSCAN2,ZFP29,ZNF854,ZSCAN2 Rabbit pAb,Zinc finger protein 29 homolog,Zinc finger protein 854

## 产品验证



Western blot analysis of ZSCAN2 expressed in human lung cancer using ZSCAN2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)