

RCOR2 Rabbit pAb

货号: **AYP20761**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	58kDa
实测分子量	58kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	histone deacetylase complex
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 250-350 of human RCOR2 (NP_775858.2).
序列	HHPLRTRRRPPKGMYSPEGLTAVSGSPDLANLTLRGLDSQLISLKRQVQSMKQTNSSLRQALEGGIDPLRPPEANTKFNS RWTTDEQLLAVQAIRRYGKD

靶点信息

研究背景	May act as a component of a corepressor complex that represses transcription.
基因ID	283248
基因名	RCOR2

Swiss	Q8IZ40
别名	RCOR2

产品验证



Western blot analysis of RCOR2 expressed in Mouse brain,Rat brain using RCOR2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn