

CAPRIN2 Rabbit pAb

货号: **AYP20789**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	126kDa
实测分子量	126kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y,293T,C6
细胞定位	centrosome,cytoplasm,cytosol,mitochondrion,nucleoplasm,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 88-329 of CAPRIN2 (NP_001002259.1).
序列	QVNHSQHGESQRALSPLQSTLSSAASPSQAYETYIENGLICLKHKIRNIEKKKLKLEDYKDRLKSGEHLNPDQLEAVEKYEE VLHNLEFAKELQKTFSGLSLDLLKAQKKAQRREHMLKLEAEKKLRTLQVQYVLQNLQEHVQKDFKGGNGAVYLP SKE LDYLIKFSKLTCPERNESLSVEDQMEQSSLYFWDLLEGSEKAVVGTTYKHLKDLLSKLLNSGYFESIPVPKNAKEKEVP

靶点信息

研究背景	The protein encoded by this gene may regulate the transport of mRNA. It may play a role in the differentiation of erythroblasts. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	65981

基因名	CAPRIN2
Swiss	Q6IMN6
别名	EEG1;EEG-1;C1QDC1;RNG140;CAPRIN2

产品验证



Western blot analysis of CAPRIN2 expressed in SH-SY5Y,293T,C6 using CAPRIN2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn