

— ABLYBIO, Help Your Research



# PAXIP1 Rabbit pAb

货号: **AYP20805**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	<b>WB:</b> 1:100 - 1:500 <b>IHC:</b> 1:50 - 1:200
理论分子量	93kDa/117kDa/121kDa
实测分子量	121KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y,Mouse kidney
细胞定位	Nucleus matrix
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 301-394 of human PAXI P1 (NP_031375.3).
------	-----------------------------------------------------------------------------------------------------------------------

## 靶点信息

研究背景	This gene is a member of the paired box (PAX) gene family and encodes a nuclear protein with six BRCT (breast cancer carboxy-terminal) domains. This protein plays a critical role in maintaining genome stability, condensation of chromatin and progression through mitosis.
基因ID	22976
基因名	PAXIP1
Swiss	Q6ZW49 ( <a href="https://www.uniprot.org/uniprotkb/Q6ZW49/entry">https://www.uniprot.org/uniprotkb/Q6ZW49/entry</a> )
别名	PAXIP1,CAGF28,CAGF29,PACIP1,PAXIP1L,PTIP,TNRC2,PAXIP1 Rabbit pAb,PAX transactivation activation domain-interacting protein

## 产品验证



Western blot analysis of PAXIP1 expressed in SH-SY5Y, Mouse kidney using PAXIP1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))