

# PAXIP1 Rabbit pAb

货号: AYP20805

## 产品信息

|       |  |
|-------|--|
| 反应    | Human,Mouse,Rat  |
| 宿主    | Rabbit   |
| 克隆性   | Polyclonal   |
| 预测反应  |  |
| 应用    | WB IHC   |
| 推荐浓度  | <b>WB:</b> 1:100 - 1:500<br><b>IHC:</b> 1:50 - 1:200   |
| 理论分子量 | 93kDa/117kDa/121kDa  |
| 实测分子量 | 121KDa   |
| 形式    | Liquid   |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated   |
| 阳性对照  | SH-SY5Y,Mouse kidney   |
| 细胞定位  | Nucleus matrix   |
| 纯化    | Affinity purification  |

## 抗原信息


|      |  |
|------|--|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 301-394 of human PAXIP1 (NP_031375.3). |
| 序列   | GNILPPEVRGNLMAAGQNLQSSERSEMIATWSPAVRTLNRITNNADIQQMNRPSNVAHILQTLSPATKNLEQQVNHSQQGHTNANAVLFSQVKV                       |

## 靶点信息

|      |  |
|------|--|
| 研究背景 | This gene is a member of the paired box (PAX) gene family and encodes a nuclear protein with six BRCT (breast cancer carboxy-terminal) domains. This protein plays a critical role in maintaining genome stability, condensation of chromatin and progression through mitosis. |
| 基因ID | 22976  |

|       |  |
|-------|--|
| 基因名   | PAXIP1   |
| Swiss | Q6ZW49   |
| 别名    | PAXIP1;CAGF28;CAGF29;PACIP1;PAXIP1L;PTIP;TNRC2 |

## 产品验证



Western blot analysis of PAXIP1 expressed in SH-SY5Y, Mouse kidney using PAXIP1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)