

— ABLYBIO, Help Your Research



CD28 Rabbit pAb

货号: **AYP20810**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	25kDa
实测分子量	50kDa/60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60,Raji,Mouse spleen,Rat spleen,Rat kidney,Rat thymus
细胞定位	cell surface,cytosol,external side of plasma membrane,immunological synapse,plasma membrane
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 180-220 of human CD28 (NP_006130.1).
------	--

靶点信息

研究背景	The protein encoded by this gene is essential for T-cell proliferation and survival, cytokine production, and T-helper type-2 development. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	940
基因名	CD28
Swiss	P10747 (https://www.uniprot.org/uniprotkb/P10747/entry)
别名	CD28,Tp44,CD28 Rabbit pAb

产品验证



Western blot analysis of CD28 expressed in Mouse spleen, Rat spleen, HL-60 using CD28 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)