

ASB2 Rabbit pAb

货号: **AYP20842**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	65kDa
实测分子量	65kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T, rat spleen, rat kidney
细胞定位	cytoskeleton, cytosol, Z disc
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 448-587 of human ASB2 (NP_057234.2).
序列	MDLGCDGEPFCFSCLYGNPAPPAPQPSSRFNDAPAADKEPSVVQFCEFVSAPEVSRWAGPIIDVLLDYVGNVQLCSRLKE HIDSFEDWAVIKEAEPPRPLAHLRVRKAIGKYRIKLLDTLPLPGLRIRYLKYENTQ

靶点信息

研究背景	This gene encodes a member of the ankyrin repeat and SOCS box-containing (ASB) protein family. These proteins play a role in protein degradation by coupling suppressor of cytokine signalling (SOCS) proteins with the elongin BC complex. The encoded protein is a subunit of a multimeric E3 ubiquitin ligase complex that mediates the degradation of actin-binding proteins. This gene plays a role in retinoic acid-induced growth inhibition and differentiation of myeloid leukemia cells. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
------	---

基因ID	51676
基因名	ASB2
Swiss	Q96Q27
别名	ASB-2;ASB2

产品验证



Western blot analysis of ASB2 expressed in 293T, rat spleen, rat kidney using ASB2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn