

GALK2 Rabbit pAb

货号: **AYP20864**

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	49kDa/50kDa
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver
细胞定位	cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 109-458 of human GALK 2 (NP_002035.1).
序列	NYFLCGLKGIQEHLGSLNLTGMNCLVDGNIPSSGLSSSSALVCCAGLVTLTLVLRNLSKVELAEICAKSERYIGTEGGGM DQISFLAEEGTAKLIEFSPLRATDVKLPSGAVFVIANSCEMKAATSHFNIRVMECLRAAKLLAKYKSLQWQDKVLRLEEV QAKLGLSLEEMLLVTEDALHPEPYNPEEICRCLGISLEELRTQILSPNTQDVLIFKLYQRAKHVYSEAAARVLQFKKICEEAPENM VQLLGLMNSHMSCRDMYECSCPELDQLVDICRKFQAQGSRLTGAGWGGCTVSMVPADKLPFLANVHKAYYQRSD GSLAPEKQSLFATKPGGGALVLEA

靶点信息

研究背景	This gene encodes a highly efficient N-acetylgalactosamine (GalNAc) kinase, which has galactokinase activity when galactose is present at high concentrations. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2014]
------	--

基因ID	2585
基因名	GALK2
Swiss	Q01415
别名	GALK2;GK2

产品验证



Western blot analysis of GALK2 expressed in Mouse liver using GALK2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn