

— ABLYBIO, Help Your Research



SNTB1 Rabbit pAb

货号: **AYP20899**

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	58kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung,Mouse liver
细胞定位	cytoplasm,cytoskeleton,focal adhesion
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human SNTB1 (NP_066301.1).
------	---

靶点信息

研究背景	Dystrophin is a large, rod-like cytoskeletal protein found at the inner surface of muscle fibers. Dystrophin is missing in Duchenne Muscular Dystrophy patients and is present in reduced amounts in Becker Muscular Dystrophy patients. The protein encoded by this gene is a peripheral membrane protein found associated with dystrophin and dystrophin-related proteins. This gene is a member of the syntrophin gene family, which contains at least two other structurally-related genes. [provided by RefSeq, Jul 2008]
基因ID	6641
基因名	SNTB1
Swiss	Q13884 (https://www.uniprot.org/uniprotkb/Q13884/entry)
别名	59-DAP,A1B,BSYN2,DAPA1B,SNT2,SNT2B1,TIP-43,SNTB1,SNTB1 Rabbit pAb,59 kDa dystrophin-associated protein A1 basic component 1,Syntrophin-2,Tax interaction protein 43

产品验证



Western blot analysis of SNTB1 expressed in Mouse lung, Mouse liver using SNTB1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)