

— ABLYBIO, Help Your Research



PEX11A Rabbit pAb

货号: **AYP20935**

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	28kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis,Mouse liver
细胞定位	peroxisome
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 110-180 of human PEX11A (NP_003838.1).
------	--

靶点信息

研究背景	This gene is a member of the PEX11 family, which is composed of membrane elongation factors involved in regulation of peroxisome maintenance and proliferation. This gene product interacts with peroxisomal membrane protein 19 and may respond to outside stimuli to increase peroxisome abundance. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Oct 2012]
基因ID	8800
基因名	PEX11A
Swiss	O75192 (https://www.uniprot.org/uniprotkb/O75192/entry)
别名	PEX11-ALPHA,PMP28,hsPEX11p,PEX11A,PEX11A Rabbit pAb,28 kDa peroxisomal integral membrane protein,Peroxin-11A,Peroxisomal biogenesis factor 11A,Protein PEX11 homolog alpha,PEX11

产品验证



Western blot analysis of PEX11A expressed in Mouse testis, Mouse liver using PEX11A Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)