

HS6ST1 Rabbit pAb

货号: **AYP20948**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	48kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,U-87MG
细胞定位	plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 302-411 of human HS6S T1 (NP_004798.3).
序列	YLFERTFNLKfirPFMQYNSTRAGGVEVDEDTIRRIEELNDLDMQLYDYAKDLFQQRYYKQRLERREQRLRSREERLLHRA KEALPREDADEPGRVPTEDYMSHIEKW

靶点信息

研究背景	The protein encoded by this gene is a member of the heparan sulfate biosynthetic enzyme family. Heparan sulfate biosynthetic enzymes are key components in generating a myriad of distinct heparan sulfate fine structures that carry out multiple biological activities. This enzyme is a type II integral membrane protein and is responsible for 6-O-sulfation of heparan sulfate. This enzyme does not share significant sequence similarity with other known sulfotransferases. A pseudogene located on chromosome 1 has been found for this gene. [provided by RefSeq, Jul 2008]
------	--

基因ID	9394
基因名	HS6ST1
Swiss	O60243
别名	HH15;HS6ST;HS6ST1

产品验证



Western blot analysis of HS6ST1 expressed in U-251MG,U-87MG using HS6ST1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn