

— ABLYBIO, Help Your Research



NUP58 Rabbit pAb

货号: **AYP20951**

产品信息

反应	Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	58kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat liver
细胞定位	nuclear envelope,nuclear membrane,nuclear pore
纯化	Affinity purification


抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 400-500 of human NUP58 (NP_054808.1).
------	--

靶点信息

研究背景	This gene encodes a member of the nucleoporin family that shares 87% sequence identity with rat nucleoporin p58. The protein is localized to the nuclear rim and is a component of the nuclear pore complex (NPC). All molecules entering or leaving the nucleus either diffuse through or are actively transported by the NPC. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]
基因ID	9818
基因名	NUP58
Swiss	Q9BVL2 (https://www.uniprot.org/uniprotkb/Q9BVL2/entry)
别名	NUP45, NUPL1, PRO2463, NUP58, NUP58 Rabbit pAb, 58 kDa nucleoporin, Nucleoporin-like protein 1, KIAA0410

产品验证



Western blot analysis of NUP58 expressed in Rat liver using NUP58 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)