

— ABLYBIO, Help Your Research



ZBPB Rabbit pAb

货号: **AYP20964**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	
实测分子量	34kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	acrosomal membrane,acrosomal vesicle,extracellular region,nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 45-200 of human ZBPB (NP_008940.2).
------	---

靶点信息

研究背景	ZPBP is one of several proteins that are thought to participate in secondary binding between acrosome-reacted sperm and the egg-specific extracellular matrix, the zona pellucida (McLeskey et al., 1998 [PubMed 9378618]).[supplied by OMIM, Aug 2008]
基因ID	11055
基因名	ZPBP
Swiss	Q9BS86 (https://www.uniprot.org/uniprotkb/Q9BS86/entry)
别名	ZPBP1,ZPBP,ZPBP Rabbit pAb,Inner acrosomal membrane IAM38,Sp38

产品验证



Western blot analysis of ZPBP expressed in Mouse brain using ZPBP Rabbit pAb at 1:1000. Secondary anti body: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)