

# STRN4 Rabbit pAb

货号: **AYP21066**

## 产品信息

反应	Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	
实测分子量	98kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Rat lung
细胞定位	cytoplasm, dendrite, dendritic spine
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 230-320 of human STRN 4 (NP_037535.2).
序列	SGGESLLVKQIEEQIKRNAAGKDGKERLGGSVLGQIPFLQNCEDSDEDELDVQHKKQRVKLPKALVPEMEDEDE EDDSEDAINEFD

## 靶点信息

研究背景	Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein.
基因ID	29888
基因名	STRN4

Swiss	Q9NRL3
别名	PPP2R6C;ZIN;zinedin;STRN4

## 产品验证



Western blot analysis of STRN4 expressed in Rat lung using STRN4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)