

— ABLYBIO, Help Your Research



# SPATS2 Rabbit pAb

货号: **AYP21072**

## 产品信息

反应	Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	
实测分子量	60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Rat testis
细胞定位	cytoplasm, cytosol
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human SPATS2 (NP_075559.2).
------	--

## 靶点信息

研究背景	Enables RNA binding activity. Located in cytosol.
------	---

基因ID	65244
基因名	SPATS2
Swiss	Q86XZ4 ( <a href="https://www.uniprot.org/uniprotkb/Q86XZ4/entry">https://www.uniprot.org/uniprotkb/Q86XZ4/entry</a> )
别名	Nbla00526,P59SCR,SCR59,SPATA10,SPATS2,SPATS2 Rabbit pAb,Serine-rich spermatocytes and round spermatid 59 kDa protein

## 产品验证



Western blot analysis of SPATS2 expressed in Rat testis using SPATS2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)